**Bürkle GmbH**

Rheinauen 5 | D - 79415 Bad Bellingen  
Telefon +49 (0) 76 35 8 27 95-0  
Telefax +49 (0) 76 35 8 27 95-31  
info@buerkle.de | www.buerkle.de

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## General safety requirements

- ▶ Read these directions for use carefully before using.
- ▶ Keep these directions for use where they will be accessible for all users.
- ▶ Failure to comply with the directions for use shall render all warranty conditions null and void and may lead to personal injury and material damage.
- ▶ The instrument must only be used by qualified skilled employees

## Safety instructions before using

**NOTE!**

The instrument is suitable in the Zones 0, 1 and 2 and it can be used for flammable liquids of explosion groups IIA and IIB (see annex).

- ▶ Follow the manufacturer's instructions and safety instructions of the medium you are filling and the pressurised gas.
- ▶ The instrument is suitable for indoor use
- ▶ Damage to the respiratory tract and other organs due to off-gassing media. Fill only in ventilated rooms or under a suitable exhaust hood.
- ▶ When working with explosive substances, follow the technical rules for combustible liquids (in accordance with TRBF – (German) Technical Rules on Combustible Liquids) and wear appropriate protective equipment.

**⚠ WARNING!**

Possible sparking due to electrostatic discharge. Danger of explosion if combustible liquids are used.

- ▶ Connect the earthing terminal before use (see scope of delivery and accessories (Antistatic set)).
- ▶ Wear the necessary protective clothing.

- ▶ The instrument must only be used within the confines of applicable material resistance and for its designated purpose. Any other type of use will affect the warranty/guarantee.

**NOTE!**

Before using, the instrument has to be cleaned to the extent of your (validated) cleaning process.

- ▶ Measures of manual cleaning have to be carried out in the cleaning/desinfection area and under consideration of occupational health and safety.
- ▶ Pay attention to the material resistance.
- ▶ Do regularly visual controls of the adapters and hoses for tightness and integrity.

## Safety instructions during and after use

- ▶ The instrument must stand on a solid conductive surface or be hung on a suitable wall.

**⚠ WARNING!**

Danger due to unsuitable condition of the workplace.

- ▶ Safe hanging/safe standing of the instrument.
- ▶ Keep instrument away from aggressive environmental conditions.
- ▶ The instrument produces excess pressure. Bodily injury due to escaping liquid. Close the shut-off valve on the discharge tube before mounting. In some types of containers loud noises can occur.

**⚠ WARNING! (→ CONTINUATION NEXT PAGE)**

In the event of a leak compressed gas can escape from all parts of the instrument (e.g. escape valve, shut-off valves or when dismantling).

**⚠ WARNING! (CONTINUATION PREVIOUS PAGE)**

- ▶ When handling nitrogen follow the producers safety data sheet.
- ▶ Ensure adequate air ventilation.
- ▶ Wear the necessary protective clothing.
- ▶ Only use the instrument under supervision.
- ▶ The instrument is not suitable for cleaning in a dishwasher.

**Layout and function**

**⚠ WARNING!**

Danger due to escaping compressed gas.

- ▶ Only use in well ventilated inner area.
- ▶ Wear the necessary protective clothing.

**⚠ WARNING!**

Some liquid will continue to run. Bodily injury due to escaping liquid.

- ▶ Close the shut-off valve discharge hose ⑥ before mounting.

- ▶ Parts in contact with medium made of stainless steel (1.4301) and PTFE.
- ▶ The instrument works with neutral gaseous media (compressed air or nitrogen) and an inlet pressure between 3 and max. 10 bar. Ensure correct connection of the supply line (see IN-access (G 1/4 "thread) or adapter (see scope of delivery) and OUT-output on the pressure regulating unit.
- ▶ Only use the instrument on UN containers made of steel or plastic (except tinplate canisters) with a test pressure of min. 150 kPa (1.5 bar, 21.7 psi).
- ▶ The instrument is suitable for filling highly pure liquids and viscous media up to 500 mPas and up to max. 200 l and flammable liquids of group IIA and IIB (see appendix). Pay attention to the chemical Resistance!
- ▶ When filling flammable liquids always earth the instrument with earthing cable of pressure regulating unit (scope of delivery) and anti-static set (accessory).
- ▶ The instrument is suitable at a temperature range from 0 to 40 °C.

**NOTE!**

- ▶ The pressure regulating unit is only suitable for use with the supplied accessories.
- ▶ Constructional changes of the instrument or the accessories aren't allowed.

**NOTE!**

Only trained specialists are allowed to carry out maintenance and servicing!

- ▶ Under normal conditions an inspection is recommended every 6 months. During this the instrument should be externally inspected for damage and checked for function. In case of intense conditions, shorter maintenance intervals may be required.
- ▶ It is recommended to regularly deflate the escape valve to ensure the function of the escape valve.


**Scope of delivery**

1 x	Pressure regulating unit	①
1 x	Transport stands	②
1 x	Pump	③
2 x	Protective Plugs	⑥
3 x	Compressed air supply adapter (1x coupling connector NW 7.2 mm, 1x hose nozzle for internal Ø 6 mm, 1 x hose nozzle for internal Ø 8 mm)	⑦
1 x	Key	⑧
1 x	Earthing cable pressure regulating unit	⑩
1 x	Connecting hose	⑫

- ▶ Check delivery for completeness and transport damage.
- ▶ Point out any damage in the packaging to the freight forwarder/supplier immediately

**Assembly**

1. Stand pressure regulating unit ① on a solid, conductive surface or hang it on a suitable wall. The hanging on the wall is possible with or without transport stand ②. The borehole spacing of instrument and transport stand is the same.

**CAUTION!**

Instrument can fall down or slip away. The operator is responsible for a safe stand or a secure wall mounting. (Screws are not included in the accessories).

2. Guide pump ③ into container and screw the barrel screw connection ④ finger-tight and air-tight.
3. If necessary open the clamping lever ⑤ and adjust the riser pipe. Then close the clamping lever and lock the riser pipe in place. Screw discharge hose on riser pipe.
4. Remove protective plugs ⑥ from IN-input and OUT-output. (IN-input: Remove by hand strength; OUT-output: Disconnect protective plug from quick coupling. To do this, first pull the slidable ring on the coupling and then disconnect the coupling. Image disconnecting quick coupling see next page.)
5. Install supply hose of pressure regulating unit with thread G 1/4" or suitable adapter ⑦. Tighten the adapter ⑦ lightly with the key ⑧. (Adapters have a sealing ring. Sealing tape isn't necessary.)
6. Turn the regulating tap ⑨ of the pressure regulating unit on (turn to the right) until the manometer indicates a pressure between 0,2 and 0,3 bar.

**! WARNING!**

Possible sparking due to electrostatic discharge. Danger of explosion if combustible liquids are used.

- ▶ Before use, connect the plug of the earthing cable pressure regulating unit ⑩ to the connection for earthing cable ⑪ and ground the pressure regulating unit. Ground pump and barrel with antistatic-set (see accessories).
- ▶ Wear the necessary protective clothing.

7. Plug connecting hose (end without connection for quick coupling) ⑫ into OUT-output of pressure regulating unit (right).

**CAUTION!**

Bursting of the hose due to abrasion or chemical stress.

- ▶ Before each use do visual control.

8. Connect connecting hose (end with connection for quick coupling) ⑬ to the quick coupling pump ⑭.

**Safe filling**

1. Open shut-off valve pump ⑮.
2. Open shut-off valve pressure regulating unit ⑯ to

produce excess pressure.

3. Hold discharge hose into filling container.

**NOTE!**


Liquid flows out when excess pressure is built up.

- ▶ The production of the excess pressure depends on the size of the barrel and its filling level and may take a while.
- ▶ To regulate the flow use regulating tap ⑯.
- ▶ If the pressure is too high, the escape valve of the pressure regulating unit opens automatically.

4. Open shut-off valve discharge hose ⑰ and fill of liquid.
5. After filling close shut-off valve discharge hose ⑰.

**Completing filling and dismantling**

1. Close shut-off valve pump ⑮.
2. Close shut-off valve pressure regulating unit ⑯.
3. Disconnect connecting hose ⑬ and quick coupling pump ⑭. To do this, first pull the slidable ring on the coupling and then disconnect the coupling.

4. Turn shut-off valve pump ⑮ on "Air". Let excess pressurised gas escape completely of the barrel. After pressure reduction, raise the discharge hose and allow the liquid to run back into the barrel.

**! WARNING!**

The container is under pressure. Bodily injury due to escaping liquid.

5. Screw off discharge hose and unscrew pump carefully from the barrel.
6. Rinse the pump with a pH-neutral fluid. If necessary, wipe the casing pressure regulating unit with a damp cleaning cloth, don't use any sharp cleaners.
7. Connect protective plugs on IN-input and OUT-output of the pressure regulating unit.

## Accessories

Type	Specification	Number	Item No.
Antistatic set	to prevent electrostatic charging	1	5602-1000

More accessories can be found at [www.buerkle.de](http://www.buerkle.de)

## Annex

Operational directive in accordance with 1999/92/EC:

Zone 0	a potentially explosive atmosphere is assumed to be continuously present
Zone 1	a potentially explosive atmosphere is assumed to be occasionally present
Zone 2	a potentially explosive atmosphere is not anticipated (and if it does occur, then only rarely and briefly)

Explosion groups:

Group I	for use in mining
Group II	for use in other areas subject to explosion hazard
Group IIA	Maximum permitted gap > 0.9 mm
Group IIB	0.5 mm < maximum permitted gap < 0.9 mm
Group IIC	Maximum permitted gap < 0.5 mm

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