**Bürkle GmbH**

Rheinauen 5 | D - 79415 Bad Bellingen  
Telefon +49 (0) 76 35 8 27 95-0  
Telefax +49 (0) 76 35 8 27 95-31  
info@buerkle.de | www.buerkle.de

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## General safety requirements

- ▶ Read these directions for use carefully before using.
- ▶ Keep these directions for use where they will be accessible for all users.
- ▶ Failure to comply with the directions for use shall render all warranty conditions null and void and may lead to personal injury and material damage.
- ▶ The instrument must only be used by qualified skilled employees.

## Safety instructions before using

- ▶ Follow the instructions and safety instructions of the gas producer.
- ▶ CO<sub>2</sub> escapes during operation. Danger of frostbite. Wear cold-protection gloves.
- ▶ CO<sub>2</sub> escapes during operation. Danger of suffocation. Fill only in ventilated rooms or under a suitable exhaust hood.

## Layout and function

The instrument is used to produce dry ice tablets weighing up to 50 g. The instrument must only be connected to CO<sub>2</sub> gas canisters (W 21.80 x 1/14" right, DIN 477 No. 6) with a riser pipe. Use only CO<sub>2</sub> gas canisters without a pressure reducer.

Point out any damage in the packaging to the freight forwarder/supplier immediately.

## Assembly

### WARNING!

CO<sub>2</sub> escapes during operation. Danger of frostbite and suffocation.

- ▶ Wear cold-protection gloves.
- ▶ Fill only in ventilated rooms or under a suitable exhaust hood.

1. Set the CO<sub>2</sub> gas canister up vertically.
2. Screw the instrument tightly onto the CO<sub>2</sub> gas canister with the hexagonal nut. The handle should be pointing down.
3. Open the tightening strap ① and remove the handle ② (arrows).
4. Open the valve of the CO<sub>2</sub> gas canister briefly to ensure the fill level of the canister and cleanliness of the evaporator nozzle.
5. Place the handle on the instrument and close the tightening strap.

## Producing dry ice tablets

### WARNING!

CO<sub>2</sub> escapes during operation. Danger of frostbite and suffocation.

- ▶ Wear cold-protection gloves.
- ▶ Fill only in ventilated rooms or under a suitable exhaust hood.

1. Slowly open the valve of the CO<sub>2</sub> gas canister.
2. Check the fill level of the instrument.
3. As soon as the tightening strap ① is stretched and the handle ② is pressed down, the maximum fill level is reached (after about 60 s). Close the valve of the CO<sub>2</sub> gas canister.
4. Hold the handle firmly and release the tightening strap.
5. Shake the handle piece a little to loosen the dry ice tablet and then remove it.

### NOTE!

The CO<sub>2</sub> gas canister cannot be completely emptied. It is designed for emptying down to a residual amount of about 20%.

## Troubleshooting

Problem	To eliminate the error
No dry ice can be produced	<ol style="list-style-type: none"> <li>1. CO<sub>2</sub> gas canister has no riser pipe.</li> <li>2. CO<sub>2</sub> gas canister riser pipe is faulty.</li> <li>3. A gas canister with a pressure reducer is used.</li> </ol>
Dry ice can't be produced correctly	<ol style="list-style-type: none"> <li>1. CO<sub>2</sub> gas canister is too hot.</li> </ol>

## Standard values

Filling weight of CO <sub>2</sub> gas canister (kg)	Tablets (number)
6,35	22
10	39
22,6	90

Average weight of the dry ice tablets: 48 g

The standard values were measured under optimum conditions. The values may vary depending on room temperature or filling duration. Cool CO<sub>2</sub> gas canisters generally produce more dry ice tablets.

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