Bürkle GmbH

Rheinauen 5 | D - 79415 Bad Bellingen
Telefon +49 (0) 76 35 8 27 95-0
Telefax +49 (0) 76 35 8 27 95-31
info@buerkle.de | www.buerkle.de

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General safety requirements

- ▶ Read these directions for use carefully before using.
- ▶ Keep these directions for use where they will be accessible for all users.
- ▶ Failure to comply with the directions for use shall render all warranty conditions null and void and may lead to personal injury and material damage.
- ▶ The instrument must only be used by qualified skilled employees.

Safety instructions before using

- ▶ Follow the instructions and safety instructions of the gas producer.
- ▶ CO₂ escapes during operation. Danger of frostbite. Wear cold-protection gloves.
- ▶ CO₂ escapes during operation. Danger of suffocation. Fill only in ventilated rooms or under a suitable exhaust hood.

Layout and function

The instrument is used to produce dry ice tablets weighing up to 250 g. The instrument must only be connected to CO₂ gas canisters (W 21.80 x 1/14" right, DIN 477 No. 6) with a riser pipe. Use only CO₂ gas canisters without a pressure reducer.

Point out any damage in the packaging to the freight forwarder/supplier immediately.

Assembly

WARNING!

CO₂ escapes during operation. Danger of frostbite and suffocation.

- ▶ Wear cold-protection gloves.
- ▶ Fill only in ventilated rooms or under a suitable exhaust hood.

1. Set the CO₂ gas canister up vertically.
2. Screw the instrument tightly onto the CO₂ gas canister with the hexagonal nut. The handle should be pointing down.
3. Open the tightening strap ① and remove the handle ② (arrows).
4. Open the valve of the CO₂ gas canister briefly to ensure the fill level of the canister and cleanliness of the evaporator nozzle.
5. Place the handle on the instrument and close the tightening strap.

Producing dry ice tablets

WARNING!

CO₂ escapes during operation. Danger of frostbite and suffocation.

- ▶ Wear cold-protection gloves.
- ▶ Fill only in ventilated rooms or under a suitable exhaust hood.

1. Slowly open the valve of the CO₂ gas canister.
2. Check the fill level of the instrument.
3. As soon as the tightening strap ① is stretched and the handle ② is pressed down, the maximum fill level is reached (after about 60 s). Close the valve of the CO₂ gas canister.
4. Hold the handle firmly and release the tightening strap.
5. Shake the handle piece a little to loosen the dry ice tablet and then remove it.

NOTE!

The CO₂ gas canister cannot be completely emptied. It is designed for emptying down to a residual amount of about 20%.

Troubleshooting

Problem	To eliminate the error
No dry ice can be produced	<ol style="list-style-type: none"> 1. CO₂ gas canister has no riser pipe. 2. CO₂ gas canister riser pipe is faulty. 3. A gas canister with a pressure reducer is used.
Dry ice can't be produced correctly	<ol style="list-style-type: none"> 1. CO₂ gas canister is too hot.

Standard values

Filling weight of CO ₂ gas canister (kg)	Tablets (number)
6,35	4
10	8
22,6	16

Average weight of the dry ice tablets: 223 g

The standard values were measured under optimum conditions. The values may vary depending on room temperature or filling duration. Cool CO₂ gas canisters generally produce more dry ice tablets.

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